
Cloud

DE Gebrauchsanweisung	4
EN Instructions for use	15
FR Instructions d'utilisation	25
IT Istruzioni per l'uso	37
ES Instrucciones de uso	48
PT Manual de utilização	59
NL Gebruiksaanwijzing	70
SV Bruksanvisning	81
DA Brugsanvisning	92
NO Bruksanvisning	102
FI Käyttöohje	112
PL Instrukcja użytkowania	122
HU Használati utasítás	134
CS Návod k použití	145
HR Upute za uporabu	156
SL Navodila za uporabo	167
SK Návod na používanie	177
TR Kullanma talimatı	188
EL Οδηγίες χρήσης	199
RU Руководство по применению	210
JA 取扱説明書	222
ZH 使用说明书	232
ZH* 使用說明書	242
KO 사용 설명서	251

1

2

3

INFORMATION

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INFORMATION

Date of last update: 2020-06-02

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

INFORMATION

- ▶ New information regarding product safety and product recalls as well as the declaration of conformity can be obtained at oa@ottobock.com or from the manufacturer's service department (see inside or outside of back cover for addresses).
- ▶ You can request this document as a PDF file at oa@ottobock.com or from the manufacturer's service department (see inside or outside of back cover for addresses). The PDF file can also be displayed in a larger size.

These instructions for use provide important information regarding the installation and proper use of the 476C00=FFC* Cloud wheelchair seat cushion (custom size: 476C00=SA200).

All users and/or their attendants must be trained by qualified personnel in the use of the product with the aid of these instructions for use. In particular, users and/or attendants must be informed of the residual risks with the aid of the safety notices.

Please note the following:

- It may be necessary to adapt the product according to the needs of the user and the medical conditions. Adjustments may only be made by qualified personnel and in close coordination with the user.
- The manufacturer reserves the right to make technical changes to the model described in these instructions for use.

2 Product description

The wheelchair seat cushion is used as a decubitus prophylaxis. It provides optimal pressure redistribution and reduces shear forces.

The seating surface consists of cells filled with the polymer gel Floam. In the middle of the seating surface, there are round Floam cells (see fig. 1, items 4, 5, 6) filled with varying densities of Floam. These can be individually positioned depending on the posture requirements. Elongated lateral Floam cells (see fig. 1, item 2) offer added postural stability and contribute to damping.

A waterproof foam base at the bottom of the seating surface has a damping effect (see fig. 1, item 3).

The inner cover (see fig. 1, item 7) closely encompasses the Floam cells and foam base and holds them in position. The filled inner cover forms the core of the seat cushion (see fig. 2) and is protected by a liquid-repellent outer cover (see fig. 1, item 8).

The components of the wheelchair seat cushion have the following functions:

Component	Figure	Function
Floam cell ¹⁾ , round, filled with polymer gel Floam, can be positioned individually		Decubitus prophylaxis, optimal pressure redistribution and reduction of shear forces
- Soft: low fill	see fig. 1, item 5	
- Medium: medium fill	see fig. 1, item 6	
- Firm: high fill	see fig. 1, item 4	
Lateral Floam cell ¹⁾ , elongated, filled with the polymer gel Floam	see fig. 1, item 2	Lateral thigh support, postural stability
Lateral inserts	see fig. 1, item 1	Serve to position the lateral Floam cells, provide postural stability
Foam base, with waterproof cover	see fig. 1, item 3	Damping
Inner cover	see fig. 1, item 7	Serves to position the Floam cells and foam components, provides postural stability
Outer cover made of PU-coated fabric, liquid-repellent, highly elastic	see fig. 1, item 8, see fig. 3	Suitable for incontinence

¹⁾ The seat cushion can be equipped with Floam cells in the respectively required number and firmness depending on the size, posture requirements and anatomy.

3 Intended use

The safe use of the product can be ensured only if it is used as intended in accordance with the information contained in these instructions for use.

3.1 Indications for use

The product is intended **exclusively** as a wheelchair seat cushion for adults with intact skin who use a wheelchair on a temporary or permanent basis.

The optimally suited seat cushion must be selected by trained, qualified personnel depending on the size and individual sitting needs of the user. It may only be used with a wheelchair in the appropriate size.

The allowed user weight is **35 kg [77 lbs] to 250 kg [551 lbs]**. We assume no liability if the product is used outside the allowed user weight range.

3.2 Indications

- Minor to pronounced or complete restrictions of mobility
- Increased decubitus risk

3.3 Contraindications

3.3.1 Absolute Contraindications

- None known

3.3.2 Relative Contraindications

- Failure to meet physical or mental requirements

3.4 Further Usage Restrictions

The product is designed for use by **one person** only. Re-use of the product is not permitted due to hygiene considerations.

3.5 Qualification

This product may only be selected and adapted by qualified personnel. Compliance with all manufacturer specifications and all applicable legal provisions is required. Please contact the manufacturer's service department for further information (see inside or outside of rear cover for addresses).

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 General safety instructions

 WARNING
Improper handling of packaging materials Risk of suffocation due to neglect of the duty to supervise ▶ Packaging materials must be kept out of the reach of children.

 CAUTION
Altered centre of gravity resulting from use of the seat cushion Risk of tipping over or falling due to unstable settings ▶ After installing the seat cushion in the wheelchair, always check the user's sitting position and the wheelchair's centre of gravity. The sitting position and centre of gravity must be appropriate for the intended use and the user's requirements. Work with the user to find the optimal sitting position. ▶ If necessary, adjust the position of the seat and axle to improve stability. Alternatively assess the suitability of anti-tip devices.

 CAUTION
Lack of testing prior to long-term use Risk of skin reddening or pressure sore formation due to incorrectly selected seat cushion ▶ Any seat cushion should be tried out by the user under medical supervision for several hours prior to long-term use. The skin should be monitored for redness and/or other damage.

4.3 Safety instructions for use

WARNING

Improper handling of ignition sources

Risk of severe burns due to seat cushions, padding and upholstery catching on fire

- ▶ The seat cushions, padding and covers fulfil the normative requirements for flame resistance. However, they may still ignite if fire is handled improperly or negligently.
- ▶ Keep away from all ignition sources, especially lit cigarettes and open flames.

WARNING

Exceeding the service life

Serious injuries due to failure to observe the manufacturer's requirements

- ▶ Using the product beyond the specified expected service life leads to increased residual risk.
- ▶ Observe the specified service life.

WARNING

Independent modification of settings

Serious injuries to the user due to improper changes to the product

- ▶ Do not modify the settings established by the qualified personnel. Only the settings described in the section "Use" in these instructions for use may be adjusted independently.
- ▶ In case of problems with the settings, please contact the qualified personnel who adjusted your product.

CAUTION

Skin damage

Skin damage or pressure points due to overloading

- ▶ Check your skin for intactness before and during use of the product.
- ▶ Pay attention to diligent skin care and pressure redistribution during interruptions in using the product.
- ▶ If skin damage or other problems occur during use, stop using the product. Consult the qualified personnel.

CAUTION

Failure to observe care instructions

Risk of pressure sores due to incorrect skin care

- ▶ All seat cushions were developed by the manufacturer to reduce pressure on the buttocks. A seat cushion alone, however, cannot eliminate sitting pressure or possible pressure sores entirely.
- ▶ When using the seat cushion, always pay attention to correct skin care – especially through a healthy diet, hygiene and regular pressure relief.

CAUTION

Extreme temperatures

Hypothermia or burns due to contact with components, failure of components

- ▶ Do not expose the product to any extreme temperatures (e.g. direct sunlight, sauna, extreme cold).
- ▶ Do not leave the product in the immediate vicinity of heaters.

⚠ CAUTION

Use of the product during diagnostic examinations and therapeutic treatment

Impairment of the examination results or the effectiveness of treatment due to interactions of the product with devices that are used

- ▶ Make sure that examinations and treatments are carried out exclusively under the prescribed conditions.

NOTICE

Improper use

Risk of damage to the cushion surface due to user error

- ▶ Do not allow the wheelchair seat cushion to come into contact with sharp objects. This also applies to animals with sharp claws such as pet cats.
- ▶ Do not expose the wheelchair seat cushion to excessive heat or open flame.
- ▶ If the wheelchair seat cushion is expected to come into contact with liquid, such as spilt drinks or episodes of incontinence, always use it in conjunction with a liquid-repellent cover.
- ▶ Only use the official Ottobock incontinence covers for the product. The incontinence covers have been optimised for the wheelchair seat cushion to maintain the pressure redistribution capability of the seat cushion. Contact your prescribing therapist, clinician or specialist dealer to obtain a spare Ottobock cover.

NOTICE

Wear of seat cushions, padding and upholstery

Loss of functionality due to prohibited continued use

- ▶ Have the seat cushions, padding and upholstery replaced promptly in case of damage.

4.4 Nameplate

Label	Meaning
	A Additional information (optional)
	B Manufacturer's logo
	C Manufacturer's reference number
	D Manufacturer's product name
	E Global Trade Item Number (DI) ¹⁾
	F Textile care instructions with symbols according to ISO 3758
	G Textile labelling
	H Manufacturer information/address
	I Symbol for medical device
	J WARNING! Read the instructions for use before using the product. Observe important safety-related information (e.g. warnings, precautions).
	K CE marking
	L Manufacturing date ²⁾
	M LOT or serial number (PI) ³⁾ incl. barcode
	N Seam
	O Sewn in, not visible: LOT or serial number Manufacturer's reference number

¹⁾ UDI-DI to GS1 standard; UDI = Unique Device Identifier, DI = Device Identifier

²⁾ YYYY = year of manufacture; MM = month of manufacture; DD = day of manufacture

³⁾ UDI-PI to GS1 standard; UDI = Unique Device Identifier, PI = Product Identifier

4.5 Side effects

The following side effects may occur during use of the product:

- Circulatory disorders, risk of pressure sores

Contact a doctor or therapist in case of problems.

5 Delivery

5.1 Scope of delivery

INFORMATION

- ▶ To order spare parts, please contact the qualified personnel.
- ▶ When ordering spare parts, please specify the width and depth of your product. Check the sewn-in label on the cover for information on the dimensions.

The product is delivered preinstalled. The scope of delivery includes:

- Seat cushion (components see Page 16)
- Instructions for use

5.2 Storage

Store the product in a dry, enclosed room with sufficient air circulation and protection from external influences. Specific information about storage conditions: see Page 25.

6 Preparing the product for use

INFORMATION

- ▶ The procedures described in this section may only be carried out by qualified personnel.
- ▶ These instructions for use must be given to the user along with the product.

INFORMATION

- ▶ The seat cushion must be set up by qualified personnel before delivery.
- ▶ Before handing over to the user or attendant, perform sufficient seat tests under medical supervision. Have the seat cushion tested by the user as needed.

6.1 Setting Up

- 1) Insert the Floam cells (see fig. 1, items 2, 4, 5, 6) as needed into the compartments in the upper section of the inner cover (see fig. 6, see fig. 2, item 1).
- 2) Insert the foam base (see fig. 1, item 3) on the bottom of the inner cover (see fig. 5).
- 3) Insert the lateral inserts (see fig. 1, item 1) into the lateral compartments of the inner cover next to the lateral Floam cells (see fig. 4).
- 4) Pull on the outer cover (see Page 22).
- 5) Position the seat cushion in the wheelchair (see Page 21).
- 6) Lower the user down carefully on to the seat cushion and carry out a trial seat fitting.

INFORMATION: Make sure the seat cushion is centred under the user and the user is comfortable in the desired sitting position.

- 7) If applicable, adapt the arrangement of the round Floam cells (see fig. 1, items 4, 5, 6) to the needs of the user.

INFORMATION: Make a sketch or take a photograph of the arrangement of the Floam cells. Provide this to the user or attendant. This allows the user or attendant to put the cells back in the intended arrangement themselves after cleaning the seat cushion.

7 Use

7.1 Positioning

⚠ CAUTION

Damage to the product, unexpected noise or odours

Health problems for the user due to improper handling

- ▶ Stop using the product immediately in the event of defects or other hazards that can result in personal injury, or if there are noises or smells that greatly differ from the delivered condition of the product.
- ▶ Contact the qualified personnel.
- ▶ Report damage to the product that may cause health problems for the user to Ottobock.

INFORMATION

- ▶ Place the seat cushion in the middle of the seat so that there is sufficient space on all sides for the wheels and frame.
- ▶ Optimum pressure redistribution is ensured when the leg support of the wheelchair is adjusted so that the thigh of the user completely rests on the seat cushion.
- ▶ Regularly remove dust and dirt particles from the hook-and-loop strap at the bottom so that the seat cushion remains firmly in place.

Correct positioning of the seat cushion in the wheelchair must be maintained for daily use:

- The label that says "BACK" (see fig. 7) points towards the back support.
- The flat side with the foam base points downwards, the side with the Floam cells points upwards.
- The hook-and-loop strap on the underside of the outer cover is used to attach the seat cushion in the wheelchair.

7.2 Functional test

Daily: loosen the seat cushion

- 1) Do not remove the outer cover for the daily functional test.
- 2) Loosen the seat cushion. To do this, press the polymer gel in the Floam cells from the sides towards the centre of the cushion, so that as much gel as possible is in the middle of the seating surface.

Monthly: check the Floam cells

- 1) Remove the outer cover for the monthly functional test (see Page 22).
- 2) Remove the foam base on the bottom of the inner cover (see fig. 5).
- 3) Note the arrangement of the Floam cells. Use the sketch or photograph prepared for you by the qualified personnel (see Page 21).

INFORMATION: Observe the safety notices at the start of the section "Removing/putting on the inner cover" (see Page 23).

- 4) Remove the Floam cells from the inner cover (see fig. 6).
- 5) Check whether the Floam cells are functional.
- 6) If necessary, have the Floam cells replaced by qualified personnel.
- 7) Insert the Floam cells and foam base (see fig. 6, see fig. 5).
- 8) Pull on the outer cover (see Page 22).

7.3 Removing/putting on the covers

7.3.1 Removing/putting on the outer cover

Removing the outer cover

The outer cover can be removed to clean the seat cushion:

- 1) Open the zip on the outer cover.
- 2) Remove the outer cover.

INFORMATION: Cleaning information see Page 23.

Putting on the outer cover

- 1) Position the cushion core with the Floam cells facing up (see fig. 2, item 1) and the foam base facing down (see fig. 2, item 2).

- 2) The BACK label on the inner cover and the BACK label on the outer cover must face in the same direction (see fig. 7).
- 3) Pull the outer cover over the cushion core with the zip on the back and the side with the hook-and-loop strap facing down.
- 4) Close the zip on the outer cover.
- 5) Ensure the outer cover fits well and there are no large folds or wrinkles on the top surface.

7.3.2 Removing/putting on the inner cover

CAUTION

Incorrectly arranged Floam cells

Pressure points due to incorrectly arranged Floam cells

- ▶ When removing the inner cover, note where you remove each Floam cell from the cover.
- ▶ When pulling on the inner cover, be sure to replace each Floam cell in the same location where it was previously removed.

Removing the inner cover

The inner cover can be removed to clean the seat cushion:

- 1) Take the lateral inserts out of the lateral compartments of the inner cover (see fig. 4).
- 2) Remove the foam base on the bottom of the cover (see fig. 5).
- 3) Note the arrangement of the Floam cells. Use the sketch or photograph prepared for you by the qualified personnel (see Page 21).

INFORMATION: In doing so, observe the safety notices at the start of the section (see Page 23).

- 4) Remove the Floam cells from the cover (see fig. 6).

INFORMATION: Cleaning information see Page 23.

Putting on the inner cover

- 1) Insert the Floam cells into the cover (see fig. 6). Replace each Floam cell in the **same location** where it was **previously removed**.
- 2) Insert the foam base on the bottom of the cover (see fig. 5).
- 3) Insert the lateral inserts into the lateral compartments of the inner cover next to the lateral Floam cells (see fig. 4).
- 4) Ensure the cover fits well and there are no large folds or wrinkles on the top surface.

7.4 Cleaning

NOTICE

Improper cleaning of the wheelchair seat cushion

Risk of damage to seat cushion due to user error

- ▶ Do not wash the lateral inserts, Floam cells or foam base in the washing machine.
- ▶ Do not dry the seat cushion and its components in the tumble dryer.
- ▶ Gently wipe the surface of the lateral inserts, Floam cells and foam base with a soft, damp cloth or sponge; do not rub. Do not use a brush or similar on the surface.
- ▶ Do not wash the covers at a temperature higher than specified. The covers can shrink when washed at higher temperatures. This could impede the seat cushion's ability to redistribute pressure and damage the fabric.
- ▶ Do not use any bleaching agents or fabric softeners for cleaning.

INFORMATION

- ▶ In case of advanced wear and tear, have the components that are available as spare parts replaced by qualified personnel.
- ▶ To disinfect the cover, add a standard hygiene detergent during washing.

- | |
|---|
| <p>▶ Use an additional incontinence cover in case of heavy incontinence. Observe the care and cleaning instructions for the incontinence cover.</p> |
|---|

Clean the product every month.

7.4.1 Cleaning the covers

- 1) Unzip and remove the outer cover (see Page 22).
- 2) Close the zip on the cover before washing.
- 3) Take the lateral inserts, foam base and Floam cells out of the inner cover (see fig. 4, see fig. 5, see fig. 6).
- 4) The covers can be washed at +60 °C [+140 °F] using a mild, environmentally friendly detergent.

Recommendation: Wash the covers using the delicate cycle at **+40 °C [+104 °F]** to prevent excessive wear and tear.

- 5) Allow to air dry.

7.4.2 Cleaning the lateral inserts, Floam cells and foam base

- 1) Take the lateral inserts, Floam cells and foam base out of the inner cover (see fig. 4, see fig. 5, see fig. 6, see Page 23).
- 2) Use lukewarm water for cleaning. If necessary, an alcohol-free disinfectant or mild detergent can be added.
- 3) Wipe the lateral inserts, Floam cells and foam base with a damp cloth or sponge.
- 4) Rinse thoroughly with clear water.
- 5) Dry the lateral inserts, Floam cells and foam base with a towel, dab or let air dry at room temperature.
- 6) Avoid exposure to direct heat (e.g. sun or heat from radiators and stoves).

8 Maintenance and repair

- The function of the product should be checked **before each use**.
- Check the product **monthly** for wear and tear and damage. The interval between inspections should be further reduced for very heavy users.
 - Check the outer and inner covers for rips, holes and other damage.
 - Verify that the zip is working properly.
 - Remove dust and dirt particles from the outer cover's hook-and-loop strap.
 - Check the lateral inserts, Floam cells and foam base for damage.
- The user is not permitted to make repairs to the product.
- Have damaged or worn parts replaced by qualified personnel as soon as possible.
- Servicing and repairs may only be carried out by qualified personnel or the manufacturer. This will ensure that only Ottobock spare parts are used for repairs.

9 Disposal information

Return the product to the qualified personnel for disposal.

All components of the product must be disposed of properly in accordance with the respective national environmental regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 Warranty

Further information regarding the warranty terms and conditions may be obtained from the qualified personnel who fitted this product, or from the manufacturer's customer service department (see inside back cover of the instructions for use for addresses).

10.3 Service life

Expected lifetime: **2 years**

The design, manufacturing and requirements for the intended use of the product are based on the expected lifetime. These also include the requirements for maintenance, ensuring effectiveness and the safety of the product.

11 Technical data

Dimensions, weights, temperatures	
Seat cushion	
Dimensions (depth x width)	See label on the cover
Weight ¹⁾ [kg (lbs)]	3.2 (7)
User weight [kg (lbs)]	35 (77)–250 (551)
Functional heights¹⁾	
Ischial zone [mm (")]	Adaptable
Support for trochanters [mm (")]	Adaptable
Thigh trough [mm (")]	Adaptable
Edge zone [mm (")]	110 (4.3)
Floam cells	
Soft (low fill) [kg (lbs)]	0.075 (0.16)
Medium (medium fill) [kg (lbs)]	0.10 (0.22)
Firm (high fill) [kg (lbs)]	0.125 (0.27)

Component	Material
Foam base	High-density PU foam
Floam (polymer gel)	Mineral oil, Kraton, acrylic
Cover of the Floam cells and foam base	PU film, double layer
Inner cover	Elastane, washable
Outer cover	PU-coated fabric, washable

¹⁾ Reference size 420 x 405 mm (16.5" x 16"), in ready-to-use condition, without incontinence cover

Ambient conditions

Temperatures and relative humidity	
Temperature during use [°C (°F)]	0 to +40 (+32 to +104)
Transport and storage temperature [°C (°F)]	0 to +40 (+32 to +104)
Relative humidity [%]	10 to 85; non-condensing

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India Pvt. Ltd.
Otto Bock Korea HealthCare Inc.
South East Asia Co., Ltd.
Otto Bock Mobility Solutions GmbH